

Date... 9/9/19

Expt. No. 6

Page No.

Monochromators in sophisticated instrument

Apparatus Available:

Laser source, Grating, Scale with measurements.

SLO:

To determine the wavelength of the given laser source using transmission diffraction grating method.

Theory -

D = the distance from the grating to the screen.

d = the spacing b/w every two lines

(same thing as every 2 sources)

If there are N lines per mm of the grating, then the space b/w every 2 adjacent lines is $d = \frac{1}{N}$

The diffraction grating formula for the principal maxima is $d \sin \theta = n \lambda$,
where n is the order of diffraction and θ = angle of diffraction.


$$\lambda = \frac{d \sin \theta}{Nn} \text{ (metre)}$$

Teacher's Signature :

Date.....

Page No.....

Expt. No.....



Date.....

Expt. No.....

Page No.....

n	$S(\text{cm})$	$2L(\text{cm})$	$l(\text{cm})$	$\tan \theta = L/s$	$\theta = \tan^{-1}(L/s)$	$\sin \theta$	Mean	$\lambda(\text{Å})$
1	25	3.2	1.6	0.064	3.66	0.0638	0.0624	637
	30	3.8	1.9	0.063	3.6	0.0627		
	35	4.5	2.25	0.064	3.66	0.0638		
2	25	6.5	3.25	0.13	7.4	0.128	0.128	640
	30	7.8	3.9	0.13	7.4	0.128		
	35	9.2	4.6	0.13	7.4	0.128		
3	25	8	4	0.16	9.09	0.157	0.207	690
	30	12	6	0.2	11.3	0.195		
	35	14.8	9.9	0.282	15.7	0.27		

Date.....

Page No.....

Expt. No.....


Observation:No. of lines per meter on the grating is 10^5 nm Calculations:

$$\frac{0.0678 + 0.0627 + 0.0678}{3} = 0.0674 \quad \lambda = \frac{0.674}{10^5} = 674$$

$$\frac{0.128 + 0.128 + 0.128}{3} = 0.128 \quad \lambda = \frac{0.128}{10^5 \times 2} = 640$$

$$\frac{0.157 + 0.195 + 0.27}{3} = 0.207 \quad \lambda = \frac{0.207}{10^5 \times 3} = 690$$

$$\frac{624 + 640 + 690}{3} = 654.6 \text{ nm}$$

Result:The wavelength of the laser source is found to be 654.6 nm

19BCE0811

9/9/19

Teacher's Signature :